Geratherm® Non-Mercury Basal Thermometer

Thank you for purchasing our "Non-Mercury" basal thermometer. The Geratherm's shell and capillary are made of medical grade glass. The liquid alloy inside (Galinstan®) is nontoxic, so if broken, clean up can be done with warm soapy water. The Geratherm non-mercury basal thermometer is a safe and effective alternative and does not require batteries that will need to be replaced.

Galinstan was developed by Geratherm Medical AG, an ISO 9002 firm located in Geschwenda, Germany, with 125 years experience of manufacturing glass thermometers.

What is a basal body thermometer (BBT)?

A basal body thermometer is an ultra sensitive thermometer that tracks your body's exact temperature.

How does this thermometer work and how accurate is the Geratherm® basal thermometer?

Just like a mercury thermometer, the Geratherm® basal thermometer measures the slightest change in body temperature. The Geratherm thermometer exceeds ASTM standards with an accuracy of +/- 0.1° F. A woman's temperature can raise an average of 0.5° F just after ovulation. This increase is caused by hormonal changes related to ovulation.

Who should use a basal thermometer?

Couples who choose to incorporate fertility awareness into their lives whether they wish to become pregnant or avoid pregnancy.

When and where should I check my temperature?

Take your temperature each morning **before you get out of bed at the same time each morning**. Place the thermometer under your tongue for 3 minutes. Don't eat, drink or smoke before you take your temperature. You may also check your temperature rectally, but choose one method to provide consistent results. Then record your temperature every day on the graph provided (back).

How do I record my temperature and why is that important?

At the top of the graph are the days of your cycle from 1-42, but add more days if needed. Underneath each cycle day you write the month and actual calendar date. Vertically listed are temperatures from 99.5 - 96.0° F. Record your temperature by using a dot and a down-pointing arrow to indicate the days you had intercourse. After ovulation, your temperature rises about 0.5° F until the beginning of the next cycle. After completing about three charts, you should notice a pattern of ovulation: temperature does not rise until after ovulation. The purpose of the chart is to assist you in observing your own individual fertility pattern.

How to begin?

The first day of your menstrual flow is day 1 of your BBT chart. Your temperature should drop when your menstrual flow starts, do not include spotting prior to your period as your first day. Record your temperature throughout your period. Make certain you note the actual day of the month in the space provided on your chart. Note any premenstrual symptoms and any true illnesses or fever and change charts when you get your period again.

What is ovulation?

Ovulation occurs when the egg is released from the ovary. Certain changes in the body may occur, such as the presence of a thin clear vaginal mucus that comes from the cervix just before ovulation, around the 10th to 12th day of your cycle. This fluid allows sperm to remain alive for days until ovulation. Other symptoms may include a slight pain on the right or left side at the time of actual ovulation. The most likely days to become pregnant, with a 28 day cycle, are the 10th and 14th day. Knowing when you are ovulating can be helpful in determining when to avoid sexual intercourse or increase the likelihood of becoming pregnant.

Physician Consultations?

It is strongly recommended that you consult with your personal physician when beginning family planning procedures. They can further assist you with interpreting the information on your tracking chart.

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