

SolarTherm[®] Speed Battery-Free Digital Thermometer Pointers



1. The **green arrow** is pointing at the On/Off button.
2. The **red arrow** shows the solar panels which charge the thermometer.
3. The **blue arrow** is pointing at the mylar sticker which is placed on top of the display and should be removed prior to use. With your fingernail, lift the sticker at the corner where it reads "LIFT" and peel it away from the display and discard.

NOTE:

- For first time use, the unit must be exposed to a light source for at least 3 hours (the solar panels should face the light source during charging and recharging).
- Press the On/Off button, you will hear a quick beep then the thermometer will do an internal check to make sure it is working properly; you will see 188.8^{°F}_M, then L^{°F}_M, then 188.8^{°F}_M again, and finally L^{°F} (with the ^{°F} blinking) indicating it is ready to take the temperature. Once the reading is complete, you will hear a series of beeps and see the reading on the display.
- When the low power symbol "☹" appears, the solar cell needs to be recharged at least 60 minutes. Turn the power off before recharging.
- The thermometer should not be stored in a drawer, closet or purse if you plan on taking your temperature within 3 hours.